
Tiramed
Home Care

INFRARED EAR AND FOREHEAD THERMOMETER

PG-IRT1603

Please read the product manual carefully before use.

USER MANUAL

Safety caution items.....2
Product introduction.....3
Use caution items.....4
Common sense about body temperature.....5
Product layout.....6
The screen display marks description.....6
Battery installation explanation.....7
Basic parameter instructions8
Introduction of Measurement methods.....8
Product cleaning methods instructions.....10
The instructions for using disposable earmuffs.....11

The warning signs and illustrations indicated in the manual are intended to enable you to use the product safely and correctly and to prevent any harm to you and others. Warning signs, illustrations and their meanings are as follows.

Legends	
	This mark means warning things(the things that you have to follow), the left graph shows general compulsion
	It means general warning.
	This mark means prohibited things (not allowed things), it means general prohibition in left picture.
	This mark means prohibited disassembly.
	Type BF Applied part
	Water resistance, moisture resistance.
	It means that you need to read the manual carefully before use.
	It means the package of this unit can comply with the requirements of green environmental protection.
	It means that the material of product or product itself is made of renewable material so that we can cycle this product, which benefits the environment and our earth.
	No discard it at will is allowed.
Ip22 : IP Classification	

Warning	
	Please ask professional doctors to explain the measured value of body temperature. The product is only used for human body temperature measurement. Please do not use this product for purposes other than body temperature measurement. No mobile phones are allowed to use around this product. Please do not use equipment that generates electromagnetic fields near the product.
	Please do not disassemble or repair this device by yourself including changing device. Please do not bend or stretch the device hard. Please do not impact or drop this product.

Warning	
	This unit is only used for human body temperature measurement without use of disease diagnosis; it cannot be used for emergency and continuous measurement in surgery. The kids should get away from this product. If the kids are in this range from zero to twelve years old, they need to use this unit to measure under their parents help. diagnose the disease and get treatment by themselves on the base of measurement result, they must follow the instruction of doctors. Children under 12 years of age and those who cannot express their thoughts are prohibited from using them. Please do not use this product for these people who suffer from otitis external, tympanitis and other ear diseases.

If you use or store this product beyond the range of specified temperature and humidity, may–be it can not reach original performance specification.
Use environment: temperature: from +10°C – +40 , humidity: from 15%RH – 93%RH°C
Storage environment: temperature: from –25°C – +55°C , humidity: from 0%RH – 93%RH

Product introduction

Intended use: Infrared Forehead Thermometer intended to measure human body temper ture by measuring ear canal or forehead.
Scope of application: It is suitable for displaying the body temperature of the measured object by measuring the heat radiation in the ear canal or forehead.

Features:
1.Non–touching type infrared measurement of ear/ forehead temperature.
2. Multiple colors and backlight display: White, Green, Orange and Red.
3. 9 sets of memory values.
4. The switch of degree Fahrenheit and degree Celsius°C.(original setting is degree Celsius°C)
5. Instant measurement within 1 second.
6. The device is supplied with a n earmuff for proper use. Detailed instructions are provided in this manual.
7. It has the function of sound on/off.
8. The machine idle time of 30 seconds, turn off power automatically

Warm tip:The measurement results of this unit at any time only can be a reference, it can not replace the medical diagnosis of professional doctor. If you have any question about the per–sonal temperature measurement result , please use and diagnose under the right instructions of doctor.

1. It is very dangerous for patients to judge and treat on their own only by measurement results, so please be sure to follow the doctor's instructions.
Self–judgment may lead to a worsening condition of patient.
2. Please do not touch with your hands or blow infrared sensor with your mouth. When the infrared sensor is damaged or dirty, it may cause abnormal measurement results.
3. If there is a temperature difference between the storage site and the measurement site, please place it at room temperature (measurement site) for about 30 minutes before next measurement. May result in incorrect measurement results.
4. Please keep this product out of the reach of children. When the child sticks to measure by himself, his ears may get hurt. If he swallows the battery or the transparent cover accidentally , please contact your doctor immediately.
5. When measuring body temperature, please do not get close to air conditioning. Avoid affecting the measurement accuracy.
6. Before and after use every time, use a cotton swab dipped 75% alcohol to wipe the probe surface. (If you see stains, fog or water on the infrared sensor glass, please use cotton swab dipped 75% alcohol for wiping the infrared sensor glass gently). If you wipe it with toilet paper or facial tissue, it will scratch the infrared sensor resulting in incorrect measurement result. Avoid cross–infection of ear disease and affect the accuracy of the measurement.
7. The product suffers mechanical damage. There is a possibility that the measurement is not right.
8. The product touches water or immerses water accidentally, please fully dry before use, especially the water on the surface of the sensor should be clean by using cotton swab. Our aim is to avoid causing safety accidents and affecting measurement accuracy.

Caution:
1. Please do not use this product for these people who suffer from otitis external, tympanitis and other ear diseases. It is possible to worsen the affected area.
2. Please do not use this product after swimming or bath or wet ears. It is possible to have low measured temperature value.
3. Do not place waste battery to danger zone. The battery may break.
4. When measuring the ear temperature of human body, the product must be operated in the ear temperature mode. Result in inaccurate measurement result.

Suggestions
1. When you tell the doctor measured body temperature value, please tell the doctor that you measured it with ear thermometer.
2. Please don't force to impact, fall, trample and shake this product.
3. Please do not disassemble, repair and modify this product.

4. Please do not allow liquid(such as alcohol, water– drop, hot water and so on) to enter the product body because of this product without water resistance.
5. The product must be kept clean in a dry place.
6. If you find any problems, please contact the sales , you can not repair the product by yourself.
7. Please do not use it under the environment of electromagnetic interference.
8. Please deal with the waste and residue of this produc at the end of the service life according to local laws and regulations.


Common sense about body temperature

The comparison of different measurement methods.
The measured values are different if we use different measurement methods. The WHO provides normal human body temperature reference values, please see below table about the specific temperature difference.


Measurement methods	Normal body temperature
Anal temperature	36.6°C – 38°C
Oral temperature	35.5°C – 37.5°C
Axillary temperature	34.7 °C – 37.3°C
Cochlear temperature	35.8°C – 38°C

The changes in human body temperature
Human belongs to constant temperature animals , the body temperature is basically constant, but it is not totally changeless, the human body temperature is constantly changing in a day, the details as follows:

At night: Lowest Body temperature is lowest because of sleep and decreasing activity.(be–low 37°C)
In the morning: Higher From warm bed to the lower temperature room in the morning, the whole body's muscles get contractions and produce heat.
At noon: Highest After lunch, human body reaches the highest temperature and the body will adjust naturally.
Three or four o'clock in the afternoon: Lower Due to physical exertion, blood sugar decreased.
In the evening: Lowest Due to the sun down, room temperature goes down.


The screen display marks description


Battery installation flow:

1. Press battery cover, the battery cover will bounce automatically.
2. Prepare two pieces of battery with the module of 1.5V AAA (number seven battery) batteries (It is recommended to use alkaline batteries), please install it into base of battery according to correct positive and negative poles.

Use your finger to press down hard on the endpiece of battery "–" for ejecting the battery.


Remove the battery cover under pressure

Insufficient power caution:

when battery voltage is lack, LCD shows"LO"sign and battery mark is always on, which means you should replace the battery

Warm prompts

If you do not use this product for a long time, please take out battery to assure its longevity. The liquid leakage of battery will harm the product; as well as pollute our environment. It is recommended to use alkaline batteries.

The way to handle scrapped battery should match the requirements of local country government and environmental protection organ.

1.Sound function: turn on/off

1) In the power–on state, press " mem " key to set the sound on or off.

2) Press "mem" button, LCD screen will show " " it means the sound works, meanwhile we will hear a short beep sound.


3) Press "mem" button again, " " will disappear, it means the sound is off.

2.The switch °F or °C

In off status , long press "mem" button for six seconds, it can switch between Fahrenheit degree (**°F**)and Celsius degree (**°C**). Wait for 8 seconds to turn off the product automatically or press " " to turn off the product directly.

3.Memory storage function

In off status, press "mem" button, the product can read and save 9 sets of measurement values in order (as below picture shows). It will turn off automatically without operating for 30 seconds or press" "button to turn off this product with your hands.

4.Back light status instructions

When measured temperature is <34°C , it shows LO with red backlight.

When measured temperature value is 34°C – 37.1 °C , it shows green backlight.


When measured temperature value is 37.2 °C– 38.1 °C , it shows orange backlight.

When measured temperature value is 38.2 °C– 43.0 °C , it shows red backlight.

When measured temperature value is >43.0 , it shows red backlight and shows HI.

Warm prompts: This function is for reference only.

Introduction of Measurement methods

Please pull back ears of your kid who is within one year old.

Ear temperature

Please pull these persons' ears back above.(over one year old kids and adults)


Ear temperature

Forehead

Forehead temperature


1.Ear temperature measurement

1–1. Remove the head cover of the forehead ther–mometer which can measure ear temperature. After attaching the earmuff, please press "⏻" button to turn on this product to enter the ear temperature measurement mode. LCD will display measurement value last time. Measurement may be performed when "👂" is not flashing. If "👂" is blinking, it indicates that the earmuff is not properly fitted or requires refitting. Then put the thermometer in your ear, the ear temperature can be measured directly by pressing " 🌡 " button.

1–2. After 1 second, you can see the measurement result.


Note: If you did not hear the beep sound, which represents the temperature measurement has not yet been completed. Please do not remove the thermom–eter probe from ear canal at this time. (If you have closed the sound indication, it will have no sound indication.)

When LCD is flashing, we are waiting for starting measurement.

Measurement result


2.Forehead temperature measurement

Install head cover of forehead thermometer which can measure forehead temperature.

1) After installing head cover of forehead thermometer, press"⏻"button to turn on this product to enter the forehead temperature measurement mode, LCD will display meas–urement value last time. Then the forehead thermometer is aimed at the forehead and measurement distance should be 30– 50mm, the forehead temperature can be measured directly bypressing " 🌡 " button.

2) After 1 second, you can see the measurement result.

When LCD is flashing, we are waiting for starting measurement.

Measurement result

After 1 second, you can see the measurement result.


Note: If you did not hear the beep sound, which representsthe temperature measurement has not yet been completed. Please do not remove the thermometer probe from forehead at this time. (If you have closed the sound indication, it will have no sound indication.)

3. Object pattern measurement

The installation of the head cover of forehead temperature has no influence to object pattern measurement.

3–1. Long press"ⓘ"button to enter the object measurement mode for 6 seconds, and then aim the thermometer at the object, the temperature of the target object can be measured directly by pressing"🌡"button.

3–. 2 After 1 second, you can see the result of the measurement.

When LCD is flashing, we are waiting for starting measurement.

Measurement result


Product cleaning methods instructions

In order to ensure the precise measurements, it recommended to clean the device after each use. Please use cotton swab to dip 75% alcohol to wipe the thermometer probe and remove the residue and dirt. We recommend you to clean the probe at least for 5 seconds.


Then wipe the other components with a soft and dry cloth for 10 seconds at least. After clean–ing, please ensure that there is no visible dirt or residue on your thermometer, Wait 10 minutes before measuring The device can be cleaned for maximum 40000 times.

Tools: Ear cover dispenser, disposable earmuffs


1. Open the bottom cover of the plastic dispenser, paying attention to the opening direction. Align the blister pack (containing disposable ear covers) with the slot and insert it, then close the bottom cover.


2. Open the front cover and pull out a sheet of disposable ear covers, and place the sheet onto the slot.

3. Insert the ear thermometer fully into the hole, then pull it out. The side slot will catch the film sheet, pull it down to tear off one ear cover.

Tiramed
Home Care

INFRARED EAR AND FOREHEAD THERMOMETER

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