
WATERSPLASH

ORAL IRRIGATOR

5013 SMART

Please read the product manual carefully before use.

USER MANUAL

Important safeguard.....3
Configuration.....5
How to use.....6
After using the appliance.....9
Troubleshooting.....10
Disposal11

All information in this section is provided to instruct you in the correct and safe operation of this device so that you can avoid injures to yourself and others as well as damage to property.

⚠ READ ALL INSTRUCTIONS BEFORE USING

Dangers:


–Toddlers, infants, people who are unable to operate should not allowed use this appliance. Please keep it out of reach of children.
– Never alter the appliance. Also, do not try to disassemble or repair it.
– Do not place or store the product where it can fall or be pulled into tub or sink.
– Do not reach for the product that has fallen into water. Unplug immdi–ately.
– Do not use the product while bathing.

WARNINGS:

This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruc–tion concerning use of the appliance in a safe way and understand the hazards involved.
Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
The appliance is only to be used with the power supply unit provided with the appliance. The battery must be removed from the appliance before it is scrapped. The appliance must be disconnected from the supply mains when removing the battery; The battery is to be disposed of safely.
Regarding how to remove the battery, please to see “Disposal” in page 7.
Do not fill the tank with water warmer than 40 °C.

If the external flexible cable or cord of this transformer is damaged, it shall be replaced by the manufacturer or his service agent or a similar qualified person in order to avoid a Hazard.
Keep the cord away from heated surfaces.
Please do not use it if you lost the nozzle.
Do not submerge the device in water.
Do not use for any purpose other than oral cleaning.
Keep nozzle away from small children. The nozzle is small enough to be swallowed and can be become lodged in the throat.
Certain types of mouthwash may damage the
equipment by causing the housing or water reservoir to crack. Therefore, do not use mouthwash in the water reservoir.
People who experience tooth or gum pain may not be able to use this de–vice and should consult their dentist.
Do not press the nozzle too firmly against your teeth or gums. Doing so may injure your gums.

SAVE THESE INSTRUCTIONS

a. Nozzle release button
b. On/Off button
c. Mode button
d. Mode indicator
e. Handle
f. Charging station
g. Nozzle
h. Identification ring
i. Nozzle suction tube
j. Filter
k. Water tank
l. Water tank cap


Note:
There may be some liquid remaining inside the main appliance before use (tank, water suction tube, etc), this is the remaining distilled water that was used for the performance (water flow) test, so there are no health concerns.

Charging the device(Fig. 1)

1. Insert the adaptor plug into the power outlet.
2. Put the device up onto right place of the charger base .
3. The “CH” indicator blinks when the device is charging, and it will turn off after it’s fully charged.

Note:
If the “CH” indicator is blinking during or after an irrigating, that means the battery need to be recharged. An full re–charging will take about 24 hours.


Fig.1



Use the water flosser


1. Attach nozzle onto the unit. (Fig. 2)
2. Open the water tank cap to fill the tank with water which temperature not exceed 40°C, then close the cap. (Fig. 3)
3. Or, Slide the water tank down off the unit, fill the tank with water which temperature not exceed 40°C, then attach the tank onto the unit (Fig. 4)
4. Press "Mode" key to select the mode you desired. During irrigating, you can change the working mode. (Fig. 5)

Normal: for normal cleaning.
Soft:for gentle cleaning of sensitive gums.
Pulse: for cleaning and massaging the gums.



5. Hold the unit in an upright position with nozzle tip point toward your teeth.
6. Keep your mouth opened slightly so that the water can flow out.
7. Press “On/Off” key to start irrigating. (Fig. 6)


1). Directing the water jet flow at a 90–degree angle to the teeth and gums
2). Directing the water jet slowly along the teeth
3). Directing the water jet along the gum line with the nozzle tip at a 90–degree angle to the gum tissue
4). Cleaning between the teeth and around any braces, crowns or bridge–workthoroughly.
8. Press “On/Off” key and then press the nozzle release button to replace other nozzle when you need.(Fig.7)
9. Once irrigating complete, press “On/Off” key again to turn off the de–vice.

After using the appliance:(Fig.8)

1. Open the water tank cap and dispose of the water in the tank.
2. Press “On/Off” key to turn on the unit, and drain the water inside the appliance.
3. Press “On/Off” key again to turn off the unit, and close the water tank cap.
4. Wipe the device with dry cloth.

Note:

To prevent bacteria propagating, make sure there is no remaining water inside the unit. If you’re not going to use this appliance for long time, be sure wipe and dry the appliance ASAP before store it.

Maintenance:Use only water or neutral cleaning agent to clean the device. Never use caustic or abrasive agents (e.g., vinegar or lime scale removers), since they may damage the device. Do not wash using hot water of over 50°C. Keep the device away from high temperature area and avoid direct sunlight.

Main unit

1. Remove the dirt by wiping using a cloth with liquid soap and wipe off using a dry cloth.
2. Do not submerge the main unit in water to clean.
3. Dry after washing with water.

Induction charger

1. Check the charger plug and cable at regular intervals,wipe the adapter plug with a cloth.
2. Open the cover of charger base, then clean and dry charger at regular intervals.

Nozzle

1. Wash with water and wipe with soft cloth.
2. Replace the nozzle every 6 months.

Water suction tube


1. Wash with water and wipe with soft cloth.
2. Do not bend, pull or twist the water suction tube.


Tank


1. Wash with water and wipe with soft cloth.
2. Dry the tank when not using for a long period.


Problem	Possible Cause	Solution
It does not work	You've just purchased the product or it has been left unused for more than 3 months.	Charge the unit for at least 8 hours.
It only works for a few minutes even after charging	The battery has reached the end of its lift cycle.	Contact service center for support.
	Charging time was insufficient.	Charge the unit for at least 8 hours.
Water pressure is very low	Water pressure is very low	Replace the nozzle
	The nozzle is deformed	
	There is dirt stuck in the filter	Clean the filter
Water does not come out	The tank is empty	Fill up tank
	The appliance is tilted too much while in use	Use with the appliance in a vertical state


The rechargeable battery inside your appliance is easy to remove for recy–cling. At the end of the product's life and prior to disposal, please remove the rechargeable battery from the handle by following the instructions below. (Like on Fig. 9) During below operation, be sure to protect your eyes, hands, fingers, and the surface on which you work.


1. Remove 2pcs rubber feet on the bottom the handle.
2. Insert a screwdriver into the slot located at the bottom of the handle and turn counter–clockwise to release the bottom cap. (Like on Fig. 10)
3. Slide the water tank down off the device and unscrew the two screws beside of the suction tube.(Fig. 11)
4. Open the top cover and put the inner out of the unit, then unscrew all the screws of the inner.(Fig. 12)
5. Cut off the battery line and pull out the batteries from the inner.(Fig. 13)

If you have recycling questions, please contact your local waste manage–ment office.

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