MahtaTebSina

MD Medical Device CLASSE I

NON STERILE

RESPI-IN-OUT

LOT 066-25

REF 1510004001

2030-08

SINGLE PATIENT

LATEX FREE

Fabbricante / Manufacturer:
Medinet S.r.l. - Via Canova 19/a
20145 Milano - Italia
www.medinetsrl.com

IT RESPI-IN-OUT®

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EN RESPI-IN-OUT®

RESPI-IN-OUT® is a single patient medical device that allows to exercise breathing by performing a natural inspiratory and expiratory movement without interruption.

Product description
RESPI-IN-OUT® is a positive pressure and volume exerciser for respiratory gymnastics. It is designed to encourage deep inhalation and exhalation. Deep breathing helps improve ventilation in the pulmonary alveoli and can facilitate the removal of mucus from the respiratory tract.
Directions
The breathing exercises executed with the aid of the RESPI-IN-OUT® strengthen the respiratory muscles and can assist in the prevention of postoperative respiratory complications often caused by hypoventilation, such as: atelectasis, mucus stagnation, and other pulmonary problems.
Contraindications
With correct use of the device using natural breathing there are no foreseen contraindications. Consult your healthcare professional if you have any questions or concerns.
Packaging contents
A) RESPI-IN-OUT®; B) Tube with mouthpiece; C) 3 coloured connectors with grid filters, characterised by holes of different dimensions (red, yellow and green).
Warning ⚠
This device must be used strictly as indicated and described by your doctor. If dizziness occurs during use, stop and rest a few minutes before continuing. In paediatric use the child must be supervised during the exercise. As with any medical device, keep out of reach of children.
Instructions for use
Remove the device from packaging, insert the red connector into the base of the RESPI-IN-OUT® and connect the tube with the mouthpiece to the device. Sit upright to promote optimal lung expansion while keeping the device in a vertical position. In fact, if the RESPI-IN-OUT® is inclined, the force required to perform the exercise is reduced, in turn reducing the efficiency and benefits of use. Close the lips around the mouthpiece and make a normal exhalation. The rising mobile indicator will indicate the pressure in cm H_2O, with a ratio of 1:1 (for example the number 7 indicated in the column corresponds to ± 7 cm H_2O). If a forced exhalation is made, the positive pressure may increase. Substituting the red connector with the yellow connector increases the ratio from 1:1 to 1:1.5; if the green connector is used the ratio doubles from 1:1 to 1:2 (for example the number 7 with the red connector corresponds to 7 cm H_2O, while using the yellow connector corresponds to 10.5 cm H_2O; and using the green connector corresponds to around 14 cm H_2O). At the end of the exhalation the mobile indicator descends. Without removing the lips from the mouthpiece, it is possible to make an inhalation, thus completing the breathing exercise. The rising ball will indicate the volume inhaled in cc with a ratio of 1:100 (for example, the number 3 in the graduated scale corresponding with the diameter of the mobile indicator indicates 300 cc). Substituting the red connector with the yellow one, the ratio increases from 1:100 to 1:150; substituting the red connector with the green one, again on medical advice, doubles the initial ratio, going from 1:100 to 1:200 (for example, when using the yellow connector the number 7 corresponds to 1,050 cc, while a value of around 1,400 cc can be reached with the green connector). Prolonging inhalation and/or exhalation by maintaining the mobile indicator in position at the end of the breath will improve the expected results. Replace connectors only under medical direction.
Cleaning of the RESPI-IN-OUT®
The RESPI-IN-OUT® is a single-patient medical device and must only be utilised for the period prescribed by your Healthcare Professional. After use, carefully wash the mouthpiece under warm running water (around 40°C) and leave to dry. When the RESPI-IN-OUT® device is not in use it should be kept in a clean, dry place. The device does not contain latex and therefore can be disposed of in the recyclable plastics receptacle.
Notes
The RESPI-IN-OUT® is a class I medical device, no sterile and not of measure; it is a respiratory incentivator and visualiser.

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